
Designed to help reduce the effects of digitally induced eye strain and other visual stress, improving vision when performing near and intermediate visual tasks.





Actual map based on Rx: Distance: Plano ADD: 0.75

How To Order ClearPro Relax:

For young children, early presbyopes

- 1. Distance Script
- 2. Distance PD
- 3. Pupil Fitting Height

Recommended Usage:

- Young Myopes to relax accommodation which may decrease myopia
- Young Hyperopes and Emmetropes who experience eye strain and visual fatigue from reading, schoolwork, computer and other near tasks
- Patients with accommodation or convergence problems (ie. problems with their focusing system and eye teaming)
- Pre-presbyopic patients

ClearPro Relax Availability

Index	Material	Power Range
1.5	Clear, Transitions Grey/Brown, NuPolar Grey/Brown/G15	-6.00/+5.00/-6.00 cyl
1.53	Trilogy	-7.00/+6.00/-6.00 cyl
1.56	Clear, Blue Guardian, Hi-Vex, Sunwise	-7.00/+6.00/-6.00 cyl
1.58	Litrex	-7.00/+6.00/-6.00 cyl
1.59	Polycarb Gradient Grey/Brown	-15.00/+7.50/-5.00 cyl
1.61	Clear, Blue Guardian, Transitions VII	-10.00/+9.00/-6.00 cyl
1.67	Clear, Sunwise	-12.00/+9.50/-6.00 cyl

Lens Info:

Add: 0.25D - 2.00D

Fitting Cross: 2mm Above **Lens Marking:** R

Available Upgrades:

- skeye multicoat
- ARMOUR hardcoat
- SLIMLITE Process

