




Diagrammatical Lens Analysis

ALL OFFICE DESIGNS

Optilux OFFICE
PREMIUM LENS SERIES

ESSENTIAL OFFICE
PROGRESSIVE LENS SERIES

START OFFICE



PLEASE NOTE: Measurements are not to scale.

- ✓ Order as reading RX
- ✓ Power decreases from reading area
- ✓ Shift area starts below lower lid
- ✓ Pupil centre = 60% of shift power
- ✓ Full shift starts above upper lid

Reading Rx: +1.25
Shift: -1.25

OPTICARE

www.opticare.com.au