

MyoEase

Myopia Control for Children




Features

- ✓ Retinal Peripheral Defocus Principle
- ✓ Weakens contrast, slows down myopia progression
- ✓ New technology, dual effect
- ✓ Multi-pinhole micro-transparent lens defocus design
- ✓ Professionally-designed and patented

Multi-pinhole micro-transparent lens defocus design

Through digital and precise micro-engraving, the high-precision micro-transparent pinholes are evenly arranged in a ring shape to form an optical defocus area, which is integrated into the physiological curvature of the retina to slow down the growth rate of the eye axis.

Micro-transparent pinhole defocus area



Central optical zoom

MyoEase is available in these indexes:

- | | |
|------|------------------------------------|
| 1.56 | Clear, Blue Guardian |
| 1.60 | Litrex Clear, Litrex Blue Guardian |
| 1.67 | Clear, Blue Guardian |
| 1.74 | Clear, Blue Guardian |



To learn more about MyoEase and our other premium lenses, please visit our website or scan the QR code.

info@opticare.com.au | www.opticare.com.au