Your gateway to comfortable vision.

START PROGRESSIVE

Balanced Visual Fields for a Comfortable and Natural Visual Experience

Start Progressive is Opticare's recommended balanced progressive lens, ideal for all kinds of office work that requires well-balanced visual fields. It is powered by the Digital Intelligence Technology (DIT) which ensures an even distribution of weight across all visual fields, creating a seamless transition to both near and far astigmatism areas. There is minimal distortion resulting to a low unwanted astigmatism.

- Balanced Progressive Lens Design
- Optimised visual fields
- Clear and comfortable view
- Budget-friendly without compromising on quality and performance
- Fully customisable to meet the individual needs of the wearer



Contact us to know more about Start Progressive and our other premium lenses. Email: info@opticare.com.au | Phone: 1800 251 852 | Fax: 1800 789 110 www.opticare.com.au



DIT (Digital Intelligence Technology)

Precision meets innovation with extraordinary optical quality and exceptionally low unwanted astigmatism.




DESIGN SPECIFICATIONS

Available corridor length:	7 mm, 9 mm			
Fitting cross position:	0 mm			
Recommended fitting height:	19 mm - 28 mm			
Minimum fitting height:	15 mm			
Inset:	1.5 mm			

RECOMMENDED FOR

Start Progressive is ideal for professionals who require a comfortable and balanced view for distance and near.

ACTUAL CYL MAPS



Availability

Index	Clear	Blue Guardian (BG)	BG Sunmatic II	Litrex	Sunmatic II	Transitions	Polarised	NuPolar
1.5	9	O	Ø		Ø	Ø	Ø	Ø
1.56	0	O	v		S			
1.61	9	O	v	Ø	S	S	\bigcirc	v
1.67	0	O	v		O	0		Ø
1.74	0	O	v		S	Ø		



Contact us to know more about Start Progressive and our other premium lenses. Email: info@opticare.com.au | Phone: 1800 251 852 | Fax: 1800 789 110 www.opticare.com.au

